

Produce *E. coli gcvR* Promoter with Oligos

The goal is to have the workshop make this new promoter with BioBrick ends.
Will use oligo assembly method.

Desired sequence:

Need biobrick prefix end with **EcoRI sticky end**

-35

GGCATTAAAAGCAAGCAGACAGAACCGTTCTGATTGTTGTATGCATGTTTTTTTTTATGC
-10 +1 RBS
TTTCCTTAAGAACAACTCACCCCTTAAAGGAATAA

Need biobrick suffix end with **Spe I sticky end**.

5' GGCATTAAAAGCAAGCAGACAGAACCGTTCTGATTGTTGTATGCATGTTTTTTTTTAT
GCTTTTCCTTAAGAACAACTCACCCCTTAAAGGAATAA 3'

32-mer 5' -AATTCGCGGCCGCTTCTAGAGGGCATTAAAAG

56a-mer 5' -CAAGCAGACAGAACCGTTCTGATTGTTGTATGCATGTTTTTTTTTATGCTTTTCCTTA

28-mer 5' -AGAACAACCTCACCCCTTAAAGGAATAAA

56b-mer 5' -CAACAATCAGAACGGTTCTGTCTGCTTTTAAATGCCCTCTAGAAGCGGCCGCG

60-mer 5' -CTAGTTTATTCCTTTAAGGGGTGAGTTGTTCTTAAAGGAAAGCATAAAAAAACATGCATA

Analysis and Results

BioBrick Restriction Site Analysis

Sequence contains no BioBrick restriction enzymes.

Oligator Results

The 4 overlap areas have melting temperatures of 68.72°C, 64.33°C, 55.54°C, 59.94°C respectively.

Oligonucleotide Hybridization Map:

AATTCGCGGCCGCTTCTAGAGGGCATTAAAAG GGCATTAAAAGCAAGCAGACAGAACCGTTCTGATTGTTGTATGCATGTTTTTTTTTATGCTTTTCCTTAAGAACAACCTCACCCCTTAAAGGAATAA A
CGGCCGCGAAGATCTC CCGTAATTTTCGTTCTGCTTGGCAAGACTAACACATACGTACAAAAAATACGAAGGAATCTGTTGAGTGGGAATTCCTTATT TGATC

Individual Oligonucleotides: [DOWNLOAD](#)

32-mer 5' -AATTCGCGGCCGCTTCTAGAGGGCATTAAAAG

56-mer 5' -CAAGCAGACAGAACCGTTCTGATTGTTGTATGCATGTTTTTTTTTATGCTTTTCCTTA

28-mer 5' -AGAACAACCTCACCCCTTAAAGGAATAAA

56-mer 5' -CAACAATCAGAACGGTTCTGTCTGCTTTTAAATGCCCTCTAGAAGCGGCCGCG

60-mer 5' -CTAGTTTATTCCTTTAAGGGGTGAGTTGTTCTTAAAGGAAAGCATAAAAAAACATGCATA

PCR Amplify from *E. coli* colony developed by Mike Waters '10

One forward primer and three reverse primers specific for the ampicillin resistance gene present on all pSB1A2 plasmids were designed. Annealing temperature is 42° C.

The forward primer sequence chosen was: (5' CAGTTACCAATGCTTAATCA '3) box1, #65
Reverse primer 1 sequence was chosen to be (5' GGCCAGATGGTAAGCCCTCC '3) box1, #64
Reverse primer 2 sequence was chosen to be (5' GCAGGACCACTTCTGCGCTC '3) box1, #63
Reverse primer 3 sequence was chosen to be (5' CTGTAGCAATGGCAACAAGTT '3) box1, #62

The PCR of each pair of forward and reverse primer will yield PCR products of known dsDNA length: 100 bps, 200 bps, and 300 bps (Figure 10).

